

货号	产品名称	应用	反应性
16659	Anti-Synaptophysin antibody [SP11]	IHC-FoFr, IHC-P, WB	Human
16660	Anti-Estrogen Receptor alpha antibody [SP1]	Flow Cyt, ICC/IF, IHC-FoFr, IHC-P, WB	Human
16661	Anti-Progesterone Receptor antibody [SP2]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Rat
16662	Anti-ErbB 2 antibody [SP3]	IHC-P	Human
16663	Anti-Cyclin D1 antibody [SP4]	Flow Cyt, ICC, ICC/IF, IHC-Fr, IHC-P, WB	Human,Mouse,Rat
16665	Anti-p53 antibody [SP5]	Flow Cyt, ICC/IF, IHC-P	Human
16667	Anti-Ki67 antibody [SP6]	Flow Cyt, ICC/IF, IHC-FoFr, IHC-Fr, IHC-P, WB	Common marmoset,Human,Mouse,Rat
16669	Anti-CD3 antibody [SP7]	Flow Cyt, IHC-FoFr, IHC-Fr, IHC-P, WB	Human,Mouse,Pig,Rat
16670	Anti-ALK antibody [SP8]	IHC-P, WB	Human
16693	Anti-Chromogranin A antibody [SP12]	IHC-P	Human
16694	Anti-Calectrin antibody [SP13]	ICC/IF, IHC-Fr, IHC-P	Human,Mouse,Rat
16695	Anti-Placental alkaline phosphatase (PLAP) antibody [SP15]	Flow Cyt, IHC-P, WB	Human,Mouse
16697	Anti-Calcitonin antibody [SP17]	ICC/IF, IHC-P	Human
16698	Anti-CD79a antibody [SP18]	Flow Cyt, IHC-P	Human
16699	Anti-CD5 antibody [SP19]	IHC-P	Human
16700	Anti-Vimentin antibody [SP20]	Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	Human
16701	Anti-COX2 / Cyclooxygenase 2 antibody [SP21]	IHC-Fr, IHC-P	Human
16702	Anti-CD23 antibody [SP23]	Flow Cyt, IHC-P, WB	Human
21700	Anti-Ki67 antibody [SP6], prediluted	IHC-FoFr, IHC-Fr, IHC-P	Human,Mouse,Rat
21703	Anti-CD3 antibody [SP7], prediluted	Flow Cyt, IHC-FoFr, IHC-Fr, IHC-P, WB	Cynomolgus monkey,Human,Mouse
27596	Anti-Progesterone Receptor antibody [SP2], prediluted	IHC-P	Human
27597	Anti-ErbB 2 antibody [SP3], prediluted	IHC-P	Human
27598	Anti-p53 antibody [SP5], prediluted	IHC-P	Human
27607	Anti-CD79a antibody [SP18], prediluted	IHC-P	Human
27608	Anti-Vimentin antibody [SP20], prediluted	Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	Human
27620	Anti-VEGFA antibody [SP28], prediluted	IHC-P	Human
64081	Anti-Cytokeratin 5 antibody [SP27]	Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	Human,Mouse,Rat
64082	Anti-PSMA antibody [SP29]	IHC-P	Human
64085	Anti-TMEM16A antibody [SP31]	Flow Cyt, ICC/IF, IHC-P, WB	Cynomolgus monkey,Human,Mouse
64088	Anti-CD20 antibody [SP32]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
64090	Anti-Cytokeratin 20 antibody [SP33]	Flow Cyt, ICC/IF, IHC-Fr, IHC-P	Human,Rat
64091	Anti-Cytokeratin 20 antibody [SP33], prediluted	IHC-P	Human
81081	Anti-RPS6 (phospho S240 + S244) antibody [SP45]	Dot, ELISA, Flow Cyt, ICC/IF, IHC-P	Human
81083	Anti-Bax antibody [SP47]	IHC-P	Human
93620	Anti-Mesothelin antibody [SP74]	IHC-P	Human
93665	Anti-Myeloperoxidase antibody [SP72]	Flow Cyt, ICC/IF, IHC-P	Human
93689	Anti-SOX2 antibody [SP76]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
93721	Anti-Uroplakin III antibody [SP73]	IHC-P	Human
93726	Anti-Tau antibody [SP70]	Flow Cyt, IHC-P	Human
93739	Anti-ErbB 3 antibody [SP71]	Flow Cyt, IHC-P	Human
93741	Anti-Cytokeratin 18 antibody [SP69]	Flow Cyt, ICC/IF, IHC-P, WB	Human
93884	Anti-Bcl2 alpha antibody [SP66] - N-terminal	IHC-P	Human
95361	Anti-MCM2 antibody [SP85]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
95363	Anti-Glypican 3 antibody [SP86]	Flow Cyt, ICC/IF, IHC-P, WB	Human
95364	Anti-Glypican 3 antibody [SP86], prediluted	IHC-P	Human
99887	Anti-TNS4/CTEN antibody [SP83]	Flow Cyt, IHC-P	Human
99889	Anti-MSH6 antibody [SP93]	Flow Cyt, ICC/IF, IHC-P, WB	Human
99891	Anti-BOB1 antibody [SP92]	IHC-P, WB	Human
99892	Anti-FOXA1 antibody [SP88]	Flow Cyt, ICC/IF, IHC-P	Human
99963	Anti-FOXP3 antibody [SP97]	IHC-P	Human
100909	Anti-CD163 antibody [SP96] - N-terminal	IHC-P	Human
101500	Anti-CD8 antibody [SP14]	Flow Cyt, IHC-P	Human
101530	Anti-CD4 antibody [SP35], prediluted	Flow Cyt, IHC-P	Human
101531	Anti-CD44 antibody [SP37]	IHC-P, WB	Human
101532	Anti-CDX2 antibody [SP54]	Flow Cyt, ICC/IF, IHC-P	Human
101533	Anti-CD43 antibody [SP55]	Flow Cyt, ICC/IF, IHC-P, WB	Human
101539	Anti-Met [c-Met] antibody [SP59]	IHC-P	Human
101544	Anti-ErbB 3 antibody [SP71], prediluted	IHC-P	Human
101547	Anti-SOX2 antibody [SP76] - N-terminal, prediluted	IHC-P	Human
101682	Anti-Cdk2 antibody [SP80]	IHC-P, WB	Human
101683	Anti-Calpain small subunit 1 antibody [SP81]	Flow Cyt, IHC-P	Human,Mouse,Rat
101684	Anti-Calpastatin antibody [SP82]	Flow Cyt, ICC/IF, IHC-P	Human
101688	Anti-Progesterone Receptor antibody [SP42]	Flow Cyt, ICC/IF, IHC-P	Human
101690	Anti-Stathmin 1 antibody [SP49]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
101691	Anti-RPS6 (phospho S235 + S236) antibody [SP50]	Dot, Flow Cyt, ICC/IF, IHC-P	Human
105225	Anti-Androgen Receptor antibody [SP107] - N-terminal	IHC-P	Human,Mouse,Rat
105226	Anti-Carbonic Anhydrase 9/CA9 antibody [SP106]	Flow Cyt, IHC-P	Human
105227	Anti-ICOS antibody [SP98]	IHC-P	Human
105228	Anti-Ikaros antibody [SP108]	Flow Cyt, IHC-P	Human
105229	Anti-IGJ antibody [SP105]	Flow Cyt, IHC-P	Human,Mouse

105389	Anti-Nestin antibody [SP103]	Flow Cyt, ICC/IF, IHC-P, WB	Human
111004	Anti-Vinculin antibody [SP117] - C-terminal	IHC-P	Human
111101	Anti-F4/80 antibody [SP115]	Flow Cyt, IHC-Fr/I, IHC-P	Mouse
111153	Anti-Thymidylate Synthase antibody [SP112]	Flow Cyt, IHC-P	Human
115777	Anti-beta Actin antibody [SP124]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
115783	Anti-Claudin 1 antibody [SP128]	Flow Cyt, IHC-P	Human
115803	Anti-S100 beta antibody [SP127]	Flow Cyt, IHC-P	Human,Mouse,Rat
115805	Anti-VEGF Receptor 2 antibody [SP123]	IHC-P	Human
115819	Anti-iNOS antibody [SP126]	Flow Cyt, IHC-P	Human,Mouse
119695	Anti-Cytokeratin 14 antibody [SP53]	Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	Human,Mouse,Rat
119697	Anti-Cytokeratin 7 antibody [SP52]	Flow Cyt, ICC/IF, IHC-P, WB	Human
130405	Anti-Syndecan-1 antibody [SP152]	ICC/IF, IHC-P	Human
130748	Anti-alpha 1 Fetoprotein antibody [SP154]	Flow Cyt, ICC/IF, IHC-P, WB	Human
130750	Anti-Luteinizing Hormone antibody [SP132]	IHC-P	Human
130751	Anti-Villin antibody [SP145]	Flow Cyt, IHC-P, IP, WB	Human,Mouse
130752	Anti-Ku80 antibody [SP130]	Flow Cyt, ICC/IF, IHC-P, WB	Human
130759	Anti-Parathyroid Hormone antibody [SP151]	IHC-P	Human
135371	Anti-PI 3 Kinase catalytic subunit alpha/PIK3CA antibody [SP139], prediluted	IHC-P	Human
135372	Anti-CD3 antibody [SP162]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
135377	Anti-CD15 antibody [SP159]	Flow Cyt, IHC-P	Species independent
135380	Anti-Cyclin E1 antibody [SP146]	IHC-P	Human
135383	Anti-RRM1 antibody [SP167]	Flow Cyt, ICC/IF, IHC-P	Human
135384	Anti-PI 3 Kinase catalytic subunit alpha/PIK3CA antibody [SP139]	Flow Cyt, ICC/IF, IHC-P	Human
135386	Anti-CD23 antibody [SP163]	Flow Cyt, ICC/IF, IHC-P, WB	Human
150299	Anti-Glucose Transporter GLUT1 antibody [SP168]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
150300	Anti-MAP4K4/NIK antibody [SP177]	Flow Cyt, IHC-P	Human
150301	Anti-alpha smooth muscle Actin antibody [SP171]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
150304	Anti-Eph receptor A2 antibody [SP169]	Flow Cyt, IHC-P	Human
150306	Anti-IL-2 Receptor alpha antibody [SP176], prediluted	IHC-P, WB	Human
178089	Anti-CD8 antibody [SP239]	IHC-P	Human
182023	Anti-ENT1 antibody [SP120]	Flow Cyt, IHC-P	Human
183308	Anti-Bcl6 antibody [SP155] - C-terminal	Flow Cyt, IHC-P	Human
183309	Anti-MAGEA1 antibody [SP194] - C-terminal	Flow Cyt, IHC-P, WB	Human
183310	Anti-C19orf2 antibody [SP215] - C-terminal	Flow Cyt, IHC-P	Human,Rat
183311	Anti-C4d antibody [SP91] - C-terminal	IHC-P	Human
183312	Anti-VPAC1 antibody [SP234]	Flow Cyt, IHC-P	Human,Rat
183316	Anti-MCSF Receptor antibody [SP211]	Flow Cyt, IHC-P	Human,Rat
183317	Anti-Cytokeratin 10 antibody [SP99] - C-terminal	IHC-Fr, IHC-P	Human,Mouse,Rat
183318	Anti-L1 Cadherin antibody [SP183] - C-terminal	IHC-P	Human
183320	Anti-MUC4 antibody [SP241] - C-terminal	IHC-P	Human
183321	Anti-GCDFP 15 antibody [SP185] - C-terminal	IHC-P, WB	Human
183322	Anti-CD14 antibody [SP192]	Flow Cyt, IHC-P, WB	Human
183324	Anti-ATM antibody [SP224]	Flow Cyt, ICC/IF, IHC-P	Human
183325	Anti-CD7 antibody [SP94]	IHC-P	Human
183326	Anti-CD38 antibody [SP149] - C-terminal	Flow Cyt, IHC-P	Human
183327	Anti-PU.1/Spi1 antibody [SP238]	Flow Cyt, ICC/IF, IHC-P	Human
183328	Anti-LMO2 antibody [SP51] - N-terminal	Flow Cyt, IHC-P, WB	Human,Mouse
183329	Anti-Cytokeratin 4 antibody [SP210] - N-terminal	IHC-P	Human
183330	Anti-Cytokeratin 17 antibody [SP95] - C-terminal	Flow Cyt, ICC/IF, IHC-P, WB	Human
183332	Anti-ALK antibody [SP144]	IHC-P, WB	Human
183333	Anti-Liver Arginase antibody [SP156]	IHC-P	Human
183334	Anti-VPAC2 antibody [SP235]	Flow Cyt, IHC-P	Human,Mouse,Rat
183336	Anti-Cytokeratin 5 antibody [SP178]	Flow Cyt, ICC/IF, IHC-P	Human
183338	Anti-Cyclin D3/CCND3 antibody [SP207]	Flow Cyt, IHC-P, WB	Human
183339	Anti-Caldesmon/CDM antibody [SP226]	Flow Cyt, IHC-P, WB	Human
183340	Anti-Nucleophosmin antibody [SP236]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
183341	Anti-TdT antibody [SP150]	Flow Cyt, IHC-P, WB	Human
183342	Anti-Neutrophil Elastase antibody [SP203]	Flow Cyt, IHC-P	Human
183343	Anti-CD79b antibody [SP240]	Flow Cyt, ICC/IF, IHC-P, WB	Human
183344	Anti-Cytokeratin 7 antibody [SP52]	Flow Cyt, IHC-P, WB	Human
183345	Anti-CD42b antibody [SP219]	Flow Cyt, IHC-Fr, IHC-P	Human,Mouse
183346	Anti-TN-X antibody [SP230]	IHC-P	Human
183347	Anti-SPINK1/P12 antibody [SP166]	IHC-P, WB	Human
183353	Anti-LYPD3 antibody [SP204]	IHC-P	Human
183354	Anti-CD16 antibody [SP175]	Flow Cyt, IHC-P	Human
183356	Anti-Sialoadhesin/CD169 antibody [SP216]	Flow Cyt, IHC-P	Human
183357	Anti-LYPD3 antibody [SP208]	IHC-P	Human
183358	Anti-CD13 antibody [SP182]	Flow Cyt, IHC-P	Human
183359	Anti-CD35 antibody [SP197]	IHC-P, WB	Human
184856	Anti-Slit2 antibody [SP122]	Flow Cyt, IHC-P	Human,Mouse
187260	Anti-Estrogen Receptor alpha antibody [SP1] - BSA and Azide free	ICC/IF, IHC-FoFr, IHC-P, WB	Human
192842	Anti-Galectin 1 antibody [SP247] - C-terminal	Flow Cyt, IHC-P, WB	Human
192845	Anti-LYPD3 antibody [SP250] - N-terminal	Flow Cyt, IHC-P	Human
192847	Anti-CD68 antibody [SP251]	Flow Cyt, IHC-P	Human

198380	Anti-Cytokeratin 18 antibody [SP69] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
198412	Anti-TMEM16A antibody [SP31] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Cynomolgus monkey,Human,Mouse
199432	Anti-CD33 antibody [SP266]	IHC-P, WB	Human
213215	Anti-CD4 antibody [SP35]	Flow Cyt, IHC-P	Human
213522	Anti-ATM antibody [SP224] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
216602	Anti-CDKN2A/p14ARF antibody [SP271]	Flow Cyt, ICC/IF, IHC-P, WB	Human
216603	Anti-SETBP1 antibody [SP312]	IHC-P	Human
217265	Anti-VISTA antibody [SP344]	IHC-P	Human
217266	Anti-VISTA antibody [SP345]	IHC-P	Human
217267	Anti-Smad4 antibody [SP306] - C-terminal	Flow Cyt, ICC/IF, IHC-P	Human,Mouse
218848	Anti-Placental alkaline phosphatase (PLAP) antibody [SP15] - BSA and Azide free	IHC-P, WB	Human,Mouse
223498	Anti-CTAG1B antibody [SP349]	Flow Cyt, IHC-P, WB	Human
223499	Anti-EGFR (phospho Y1068) antibody [SP353]	Flow Cyt, ICC/IF, IHC-P, WB	Human
223500	Anti-ERK1 + ERK2 (phospho T202 + Y204) antibody [SP327]	Flow Cyt, IHC-P, WB	Human
223501	Anti-Gastrin antibody [SP318]	IHC-P	Human
223502	Anti-elf4G1 (phospho S1108) antibody [SP351]	IHC-P	Human
223503	Anti-Lysozyme antibody [SP329]	Flow Cyt, IHC-P, WB	Human
223504	Anti-Mamaglobin A antibody [SP319] - C-terminal	Flow Cyt, IHC-P	Human
224798	Anti-Inhibin alpha antibody [SP323]	IHC-P	Human,Mouse
224799	Anti-Prostate Specific Antigen antibody [SP317]	IHC-P, WB	Human
224800	Anti-CD11b antibody [SP331] - C-terminal	IHC-P, WB	Human
224801	Anti-Wilms Tumor Protein antibody [SP320] - N-terminal	Flow Cyt, IHC-P	Human,Mouse,Rat
224802	Anti-PAP antibody [SP313]	Flow Cyt, IHC-P	Human
224803	Anti-beta Catenin antibody [SP328]	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
224804	Anti-smooth muscle Myosin heavy chain 11 antibody [SP314] - N-terminal	Flow Cyt, IHC-Fr, IHC-P	Human,Mouse,Rat
224805	Anti-CD11b antibody [SP330] - C-terminal	Flow Cyt, IHC-P, WB	Human
224806	Anti-Wilms Tumor Protein antibody [SP321]	Flow Cyt, IHC-P, WB	Human,Mouse,Rat
225715	Anti-Cleaved PARP1 antibody [SP276]	IHC-P, WB	Human
225716	Anti-SPARC antibody [SP205]	Flow Cyt, ICC/IF, IHC-P, WB	Human
225717	Anti-CD19 antibody [SP110] - C-terminal	Flow Cyt, IHC-P	Human
225718	Anti-Tcl1 antibody [SP248]	Flow Cyt, ICC/IF, IHC-P, WB	Human
225719	Anti-N Cadherin antibody [SP90] - C-terminal	Flow Cyt, IHC-P	Human,Mouse,Rat
225720	Anti-PI 3 Kinase p85 alpha antibody [SP62]	Flow Cyt, ICC/IF, IHC-P	Human,Mouse,Rat
225721	Anti-ERCC1 antibody [SP68]	IHC-P	Human,Mouse,Rat
225722	Anti-Factor XIIIa antibody [SP196]	IHC-P	Human
226756	Anti-Sall4 antibody [SP289]	Flow Cyt, ICC/IF, IHC-P	Human
226763	Anti-Caspase-7 antibody [SP284] - C-terminal	IHC-P, WB	Human
226834	Anti-SATB1 antibody [SP287] - C-terminal	IHC-P	Human
226837	Anti-CD47 antibody [SP279]	IHC-P	Human
226839	Anti-Alpha B Crystallin antibody [SP288]	IHC-P	Human
227635	Anti-PAX5 antibody [SP34]	Flow Cyt, IHC-P	Human,Rat
227636	Anti-Caveolin-1 antibody [SP43] - C-terminal	IHC-P	Human
227637	Anti-Met (c-Met) antibody [SP44] - C-terminal	Flow Cyt, IHC-P	Human
227638	Anti-delta 1 Catenin/CAS antibody [SP63] - C-terminal	Flow Cyt, ICC/IF, IHC-P	Human
227639	Anti-E Cadherin antibody [SP64]	IHC-P, WB	Human
227640	Anti-CD10 antibody [SP67]	Flow Cyt, IHC-P, WB	Human
227641	Anti-GFAP antibody [SP78] - C-terminal	IHC-P	Human
227642	Anti-EGFR antibody [SP84] - C-terminal	Flow Cyt, ICC/IF, IHC-P, WB	Human
227643	Anti-KRT6C antibody [SP87]	Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	Human,Mouse,Rat
227644	Anti-Cytokeratin 8 antibody [SP102] - C-terminal	Flow Cyt, ICC/IF, IHC-Fr, IHC-P	Human
227645	Anti-CD22 antibody [SP104] - C-terminal	IHC-P	Human
227646	Anti-MUM1 antibody [SP114] - C-terminal	Flow Cyt, IHC-P, WB	Human
227647	Anti-EGFR (mutated L858 R) antibody [SP125]	Flow Cyt, IHC-P, WB	Human
227648	Anti-KMT6 / EZH2 antibody [SP129] - C-terminal	Flow Cyt, IHC-P, WB	Human
227649	Anti-FOXP1 antibody [SP133] - C-terminal	Flow Cyt, IHC-Fr, IHC-P, WB	Human,Mouse,Rat
227650	Anti-TWEAKR/FN14 antibody [SP134] - C-terminal	Flow Cyt, IHC-P, WB	Human
227651	Anti-Desmin antibody [SP138] - C-terminal	ICC/IF, IHC-P, WB	Human,Mouse,Rat
227652	Anti-TTF1 antibody [SP141]	Flow Cyt, IHC-P	Human,Mouse,Rat
227653	Anti-Lambda Light chain antibody [SP147]	IHC-P, WB	Human
227654	Anti-Kappa light chain antibody [SP148]	IHC-P, WB	Human
227655	Anti-p53 antibody [SP161]	Flow Cyt, ICC/IF, IHC-P	Human
227656	Anti-CD31 antibody [SP164] - C-terminal	Flow Cyt, IHC-P	Human
227657	Anti-PTEN antibody [SP170] - N-terminal	Flow Cyt, IHC-P, WB	Human
227658	Anti-NRAS (mutated Q61 R) antibody [SP174]	Flow Cyt, IHC-P, WB	Human
227659	Anti-CD10 antibody [SP179] - N-terminal	Flow Cyt, IHC-P, WB	Human,Mouse,Rat
227660	Anti-Calponin 1 antibody [SP180] - C-terminal	IHC-P	Human
227661	Anti-Calponin 1 antibody [SP181]	IHC-P, WB	Human
227662	Anti-CD21 antibody [SP186] - C-terminal	Flow Cyt, IHC-P, WB	Human
227663	Anti-CD13 antibody [SP187]	Flow Cyt, ICC/IF, IHC-P, WB	Human
227664	Anti-MAGEA1 antibody [SP188]	Flow Cyt, IHC-P, WB	Human
227665	Anti-CD16 antibody [SP189] - C-terminal	Flow Cyt, IHC-P, WB	Human
227666	Anti-CD35 antibody [SP191]	IHC-P	Human

227667	Anti-Calponin 1 antibody [SP193] - C-terminal	IHC-P	Human,Mouse,Rat
227668	Anti-CD21 antibody [SP199]	Flow Cyt, IHC-P, WB	Human
227669	Anti-CD42b antibody [SP202]	IHC-Fr, IHC-P	Human,Mouse,Rat
227670	Anti-CD276 antibody [SP206]	IHC-P, WB	Human
227671	Anti-Hemoglobin subunit alpha antibody [SP212]	IHC-P, WB	Human
227672	Anti-Sialoadhesin/CD169 antibody [SP213] - C-terminal	Flow Cyt, IHC-P	Human
227673	Anti-PHLPP2 antibody [SP220]	IHC-P, WB	Human
227674	Anti-PTEN antibody [SP227]	IHC-P, WB	Human
227675	Anti-HLA-DP antibody [SP228]	Flow Cyt, ICC/IF, IHC-P	Human
227676	Anti-HLA-DPB1 antibody [SP229]	Flow Cyt, ICC/IF, IHC-P	Human
227677	Anti-Fbxw7 antibody [SP237]	Flow Cyt, IHC-P	Human
227678	Anti-Androgen Receptor antibody [SP242] - C-terminal	IHC-P	Human
227679	Anti-CD276 antibody [SP265] - C-terminal	Flow Cyt, ICC/IF, IHC-P, WB	Human
227680	Anti-SOX10 antibody [SP267]	Flow Cyt, ICC/IF, IHC-P, WB	Human
227681	Anti-PD1 antibody [SP269]	Flow Cyt, IHC-P	Human
227682	Anti-PDK1 antibody [SP270] - C-terminal	IHC-P	Human
227683	Anti-FGFR2 antibody [SP273] - N-terminal	Flow Cyt, IHC-P	Human
227684	Anti-SOX10 antibody [SP275]	Flow Cyt, IHC-P, WB	Human
227685	Anti-SATB1 antibody [SP280] - C-terminal	Flow Cyt, IHC-P	Human
227686	Anti-SATB2 antibody [SP281] - C-terminal	IHC-FoFr, IHC-P, WB	Human,Mouse,Rat
227687	Anti-STAT5 antibody [SP283] - C-terminal	IHC-P	Human
227688	Anti-CD19 antibody [SP291] - C-terminal	Flow Cyt, IHC-P, WB	Human
227689	Anti-TROP2 antibody [SP293] - C-terminal	Flow Cyt, IHC-P, WB	Human
227690	Anti-TROP2 antibody [SP294]	IHC-P, WB	Human
227691	Anti-TROP2 antibody [SP295]	IHC-P, WB	Human
227694	Anti-ALDH3A1 antibody [SP298] - C-terminal	IHC-P, WB	Human
227695	Anti-GPNMB antibody [SP299]	Flow Cyt, IHC-P, WB	Human
227696	Anti-SLC39A6/ZIP-6 antibody [SP300]	Flow Cyt, IHC-P, WB	Human
227697	Anti-CD103 antibody [SP301]	Flow Cyt, IHC-P	Human
227698	Anti-CD2 antibody [SP304]	IHC-P	Human
227699	Anti-TRBC1 antibody [SP305]	IHC-P, WB	Human
227701	Anti-Thrombomodulin antibody [SP309] - C-terminal	Flow Cyt, IHC-P	Human
227702	Anti-Integrin beta 3 antibody [SP324]	IHC-P, WB	Human
227703	Anti-Fibroblast activation protein, alpha antibody [SP325]	IHC-P	Human
227704	Anti-TMEM173 antibody [SP338]	Flow Cyt, ICC/IF, IHC-P, WB	Human
227705	Anti-TMEM173 antibody [SP339]	Flow Cyt, IHC-P, WB	Human
227706	Anti-LAG-3 antibody [SP346] - C-terminal	IHC-P	Human
227707	Anti-PAX8 antibody [SP348] - N-terminal	Flow Cyt, ICC/IF, IHC-P, WB	Human
227708	Anti-Lysozyme antibody [SP350] - C-terminal	Flow Cyt, IHC-P, WB	Human
227709	Anti-CTLA4 antibody [SP355]	IHC-P	Human
227737	Anti-EGFR (deleted E746 + L747 + R748 + E749 + A750) antibody [SP111]	Flow Cyt, ICC/IF, IHC-P, WB	Human
227738	Anti-CD99 antibody [SP119]	IHC-P	Human
228137	Anti-Progesterone Receptor antibody [SP42], prediluted	IHC-P	Human
228139	Anti-Stathmin 1 antibody [SP49], prediluted	IHC-P	Human
228141	Anti-CDX2 antibody [SP54] - N-terminal, prediluted	IHC-P	Human
228144	Anti-Mesothelin antibody [SP74] - C-terminal, prediluted	IHC-P	Human
228149	Anti-Cytokeratin 8 antibody [SP102], prediluted	IHC-P	Human
228151	Anti-Androgen Receptor antibody [SP107], prediluted	IHC-P	Human
228154	Anti-MUM1 antibody [SP114], prediluted	IHC-P	Human
228156	Anti-EGFR (mutated L858 R) antibody [SP125], prediluted	IHC-P	Human
228161	Anti-CD38 antibody [SP149], prediluted	IHC-P	Human
228176	Anti-CD42b antibody [SP202], prediluted	IHC-P	Human
228178	Anti-CD276 antibody [SP206], prediluted	IHC-P	Human
228180	Anti-MCSF Receptor antibody [SP211], prediluted	IHC-P	Human
228184	Anti-Galectin 1 antibody [SP247] - C-terminal, prediluted	Flow Cyt, ICC/IF, IHC-P	Human
228185	Anti-CD68 antibody [SP251], prediluted	IHC-P	Human
228188	Anti-SOX10 antibody [SP267] - C-terminal, prediluted	IHC-P	Human
228193	Anti-CD47 antibody [SP279], prediluted	IHC-P	Human
228199	Anti-TROP2 antibody [SP295], prediluted	IHC-P	Human
228205	Anti-Smad4 antibody [SP306] - C-terminal, prediluted	IHC-P	Human
228215	Anti-Fibroblast activation protein, alpha antibody [SP325], prediluted	IHC-P	Human
228229	Anti-CTLA4 antibody [SP355], prediluted	IHC-P	Human
228278	Anti-CD16 antibody [SP189] - C-terminal, prediluted	IHC-P	Human
228289	Anti-LYPD3 antibody [SP208], prediluted	IHC-P	Human
228462	Anti-PD-L1 antibody [SP142] - C-terminal	ICC/IF, IHC-P	Chinese hamster,Human
228463	Anti-PD-L1 antibody [SP142] - C-terminal, prediluted	IHC-P	Human
228464	Anti-CD134 / OX40L receptor antibody [SP195] - C-terminal	IHC-P	Human
228466	Anti-PTEN antibody [SP218] - C-terminal	Flow Cyt, IHC-P, WB	Human
228467	Anti-PTEN antibody [SP218] - C-terminal, prediluted	IHC-P	Human
228468	Anti-Indoleamine 2, 3-dioxygenase antibody [SP260] - C-terminal	Flow Cyt, ICC/IF, IHC-P	Human
228470	Anti-Indoleamine 2, 3-dioxygenase antibody [SP277] - N-terminal	IHC-P	Human
228472	Anti-Indoleamine 2, 3-dioxygenase antibody [SP282]	IHC-P	Human
228474	Anti-SLC34A2 antibody [SP322] - N-terminal	Flow Cyt, ICC/IF, IHC-P	Human
230903	Anti-CD14 antibody [SP192] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human

230904	Anti-CD7 antibody [SP94] - BSA and Azide free	IHC-P	Human
231172	Anti-Ki67 antibody [SP6] - BSA free	Flow Cyt, ICC/IF, IHC-FoFr, IHC-Fr, IHC-P, WB	Common marmoset,Human,Mouse,Rat
231434	Anti-EGFR antibody [SP9]	Flow Cyt, IHC-P	Human
231436	Anti-CD31 antibody [SP38]	IHC-P	Human
231437	Anti-MSH2 antibody [SP46]	IHC-P	Human
231438	Anti-ErbB 2 antibody [SP101]	IHC-P	Human
231439	Anti-p57 Kip2 antibody [SP118]	IHC-P	Human
231440	Anti-CD1a antibody [SP157]	IHC-P	Human
231441	Anti-IL-2 Receptor alpha antibody [SP176]	IHC-P, WB	Human
236216	Anti-Cytokeratin 5 antibody [SP27] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
236217	Anti-Mesothelin antibody [SP74] - BSA and Azide free	IHC-P	Human
236218	Anti-Myeloperoxidase antibody [SP72] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
236219	Anti-Tau antibody [SP70] - BSA and Azide free	Flow Cyt, IHC-P	Human
236220	Anti-ErbB 3 antibody [SP71] - BSA and Azide free	Flow Cyt, IHC-P	Human
236221	Anti-Bcl2 alpha antibody [SP66] - BSA and Azide free	Flow Cyt, IHC-P	Human
236222	Anti-Progesterone Receptor antibody [SP42] - BSA and Azide free	IHC-P	Human
236223	Anti-Stathmin 1 antibody [SP49] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
236224	Anti-RPS6 (phospho S235 + S236) antibody [SP50] - BSA and Azide free	Dot, Flow Cyt, ICC/IF, IHC-P	Human
236225	Anti-Androgen Receptor antibody [SP107] - BSA and Azide free	IHC-P	Human
236226	Anti-ICOS antibody [SP98] - BSA and Azide free	IHC-P	Human
236227	Anti-Ikaros antibody [SP108] - BSA and Azide free	Flow Cyt, IHC-P	Human
236228	Anti-Luteinizing Hormone antibody [SP132] - BSA and Azide free	IHC-P	Human
236229	Anti-Parathyroid Hormone antibody [SP151] - BSA and Azide free	IHC-P	Human
236230	Anti-CD8 antibody [SP239] - BSA and Azide free	IHC-P	Human
236231	Anti-C4d antibody [SP91] - BSA and Azide free	IHC-P	Human
236232	Anti-CD38 antibody [SP149] - BSA and Azide free	Flow Cyt, IHC-P	Human
236233	Anti-Liver Arginase antibody [SP156] - BSA and Azide free	IHC-P	Human
236234	Anti-CD47 antibody [SP279] - BSA and Azide free	IHC-P	Human
236235	Anti-CD10 antibody [SP67] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
236236	Anti-CD35 antibody [SP191] - BSA and Azide free	IHC-P	Human
236237	Anti-Androgen Receptor antibody [SP242] - BSA and Azide free	IHC-P	Human
236238	Anti-PD-L1 antibody [SP142] - BSA and Azide free	IHC-P	Human
236239	Anti-MSH2 antibody [SP46] - BSA and Azide free	Flow Cyt, IHC-P	Human
236434	Anti-CD20 antibody [SP32] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
236435	Anti-Cytokeratin 20 antibody [SP33] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-Fr, IHC-P	Human,Rat
236436	Anti-CD44 antibody [SP37] - BSA and Azide free	IHC-P, WB	Human
236437	Anti-Caveolin-1 antibody [SP43] - BSA and Azide free	IHC-P	Human
236438	Anti-RPS6 (phospho S240 + S244) antibody [SP45] - BSA and Azide free	Dot, ELISA, Flow Cyt, ICC/IF, IHC-P	Human
236439	Anti-Cytokeratin 14 antibody [SP53] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-Fr, IHC-P, WB	Human,Mouse,Rat
236440	Anti-Pi 3 Kinase p85 alpha antibody [SP62] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human,Mouse,Rat
236441	Anti-delta 1 Catenin/CAS antibody [SP63] - BSA and Azide free	IHC-P	Human
236442	Anti-ERCC1 antibody [SP68] - BSA and Azide free	IHC-P	Human
236444	Anti-CD163 antibody [SP96] - BSA and Azide free	IHC-P	Human
237772	Anti-CD19 antibody [SP291] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
237773	Anti-PAX5 antibody [SP34] - BSA and Azide free	Flow Cyt, IHC-P	Human
238271	Anti-CTLA4 antibody [SP355] - BSA and Azide free	IHC-P	Human
238272	Anti-IL2 Receptor alpha antibody [SP176] - BSA and Azide free	IHC-P, WB	Human
238773	Anti-EGFR (mutated L858 R) antibody [SP125] - BSA and Azide free	IHC-P, WB	Human
238776	Anti-GCDFFP 15 antibody [SP185] - BSA and Azide free	IHC-P, WB	Human
238782	Anti-Caldesmon/CDM antibody [SP226] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
238784	Anti-CD33 antibody [SP266] - BSA and Azide free	IHC-P, WB	Human
238785	Anti-Indoleamine 2, 3-dioxygenase antibody [SP277] - BSA and Azide free	IHC-P	Human
238787	Anti-SATB2 antibody [SP281] - BSA and Azide free	Flow Cyt, IHC-FoFr, IHC-P, WB	Human
238788	Anti-Indoleamine 2, 3-dioxygenase antibody [SP282] - BSA and Azide free	IHC-P	Human
238789	Anti-TROP2 antibody [SP294] - BSA and Azide free	IHC-P, WB	Human
238791	Anti-TROP2 antibody [SP295] - BSA and Azide free	IHC-P, WB	Human
238793	Anti-SLC34A2 antibody [SP322] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
238794	Anti-CD11b antibody [SP331] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
238795	Anti-TMEM173 antibody [SP338] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
238796	Anti-TMEM173 antibody [SP339] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
238798	Anti-CD4 antibody [SP35] - BSA and Azide free	Flow Cyt, IHC-P	Human
238799	Anti-CD43 antibody [SP55] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
238801	Anti-Calpain small subunit 1 antibody [SP81] - BSA and Azide free	Flow Cyt, IHC-P	Human,Mouse,Rat
238802	Anti-TNS4/CTEN antibody [SP83] - BSA and Azide free	Flow Cyt, IHC-P	Human
238803	Anti-EGFR antibody [SP84] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
238804	Anti-Glypican 3 antibody [SP86] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
238805	Anti-MSH6 antibody [SP93] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
238808	Anti-Cytokeratin 17 antibody [SP95] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
238809	Anti-FOXP3 antibody [SP97] - BSA and Azide free	IHC-P	Human
238810	Anti-Cytokeratin 10 antibody [SP99] - BSA and Azide free	IHC-P	Human
239751	Anti-SOX10 antibody [SP275] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
239793	Anti-Progesterone Receptor antibody [SP2] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human,Rat
239794	Anti-Cyclin D1 antibody [SP4] - BSA and Azide free	Flow Cyt, ICC, ICC/IF, IHC-P, WB	Human,Mouse,Rat

239795	Anti-Calretinin antibody [SP13] - BSA and Azide free	ICC/IF, IHC-Fr, IHC-P	Human,Mouse,Rat
239796	Anti-CD5 antibody [SP19] - BSA and Azide free	IHC-P	Human
239921	Anti-BOB1 antibody [SP92] - BSA and Azide free	IHC-P, WB	Human
239990	Anti-iNOS antibody [SP126] - BSA and Azide free	IHC-P	Human,Mouse
240092	Anti-CD15 antibody [SP159] - BSA and Azide free	Flow Cyt, IHC-P	Species independent
240093	Anti-RRM1 antibody [SP167] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
240258	Anti-ENT1 antibody [SP120] - BSA and Azide free	Flow Cyt, IHC-P	Human
240265	Anti-MCSF Receptor antibody [SP211] - BSA and Azide free	Flow Cyt, IHC-P	Human
240266	Anti-L1 Cadherin antibody [SP183] - BSA and Azide free	IHC-P	Human
240267	Anti-LMO2 antibody [SP51] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human,Mouse
240268	Anti-CD42b antibody [SP219] - BSA and Azide free	Flow Cyt, IHC-Fr, IHC-P	Human,Mouse
240402	Anti-CD19 antibody [SP110] - BSA and Azide free	Flow Cyt, IHC-P	Human
240403	Anti-N Cadherin antibody [SP90] - BSA and Azide free	Flow Cyt, IHC-P	Human
240404	Anti-CD22 antibody [SP104] - BSA and Azide free	Flow Cyt, IHC-P	Human
240405	Anti-Kappa light chain antibody [SP148] - BSA and Azide free	IHC-P, WB	Human
240406	Anti-CD42b antibody [SP202] - BSA and Azide free	IHC-P	Human
240407	Anti-CD276 antibody [SP206] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
240408	Anti-Hemoglobin subunit alpha antibody [SP212] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
240409	Anti-EGFR (deleted E746 + L747 + R748 + E749 + A750) antibody [SP111] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
240410	Anti-ErbB 2 antibody [SP101] - BSA and Azide free	IHC-P	Human
240411	Anti-CD1a antibody [SP157] - BSA and Azide free	IHC-P	Human
240675	Anti-alpha 1 Fetoprotein antibody [SP154] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
240927	Anti-Syntaxin 1a antibody [SP36] - BSA and Azide free	Flow Cyt, IHC-P	Human
240930	Anti-Bax antibody [SP47] - BSA and Azide free	Flow Cyt, IHC-P	Human
240932	Anti-Uroplakin III antibody [SP73] - BSA and Azide free	IHC-P	Human
240933	Anti-MCM2 antibody [SP85] - BSA and Azide free	Flow Cyt, IHC-P, WB	Dog,Human,Mouse,Rat
240935	Anti-FOXA1 antibody [SP88] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
240936	Anti-CDX2 antibody [SP54] - BSA and Azide free	Flow Cyt, IHC-P	Human
240937	Anti-Cdk2 antibody [SP80] - BSA and Azide free	IHC-P, WB	Human
240938	Anti-Calpastatin antibody [SP82] - BSA and Azide free	Flow Cyt, IHC-P	Human
240939	Anti-IGJ antibody [SP105] - BSA and Azide free	Flow Cyt, IHC-P	Human
240945	Anti-Vinculin antibody [SP117] - BSA and Azide free	IHC-P	Human
240946	Anti-F4/80 antibody [SP115] - BSA and Azide free	Flow Cyt, IHC-Fr, IHC-P	Mouse
240947	Anti-Thymidylate Synthase antibody [SP112] - BSA and Azide free	IHC-P	Human
240961	Anti-CD35 antibody [SP197] - BSA and Azide free	IHC-P, WB	Human
240972	Anti-Galectin 1 antibody [SP247] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
240980	Anti-Gastrin antibody [SP318] - BSA and Azide free	IHC-P	Human
240981	Anti-Mammaglobin A antibody [SP319] - BSA and Azide free	Flow Cyt, IHC-P	Human
240982	Anti-Prostate Specific Antigen antibody [SP317] - BSA and Azide free	IHC-P, WB	Human
240983	Anti-smooth muscle Myosin heavy chain 11 antibody [SP314] - BSA and Azide free	Flow Cyt, IHC-Fr, IHC-P	Human,Mouse,Rat
240984	Anti-E Cadherin antibody [SP64] - BSA and Azide free	IHC-P, WB	Human
240985	Anti-KRT6C antibody [SP87] - BSA and Azide free	Flow Cyt, IHC-P, WB	Cow,Human
240986	Anti-Cytokeratin 8 antibody [SP102] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-Fr, IHC-P	Human
240987	Anti-CD21 antibody [SP186] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
240988	Anti-CD21 antibody [SP199] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
240989	Anti-Fc gamma R1 antibody [SP320] - BSA and Azide free	IHC-P	Human
240990	Anti-CD99 antibody [SP119] - BSA and Azide free	IHC-P	Human
241407	Anti-LAG-3 antibody [SP346] - BSA and Azide free	IHC-P	Human
241408	Anti-CD11b antibody [SP330] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
242379	Anti-p53 antibody [SP5] - BSA and Azide free	IHC-P	Human
242380	Anti-CD23 antibody [SP23] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
242387	Anti-beta Actin antibody [SP124] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
242388	Anti-S100 beta antibody [SP127] - BSA and Azide free	Flow Cyt, IHC-P	Human
242394	Anti-Syndecan-1 antibody [SP152] - BSA and Azide free	ICC/IF, IHC-P	Human
242395	Anti-alpha smooth muscle Actin antibody [SP171] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
242396	Anti-Eph receptor A2 antibody [SP169] - BSA and Azide free	Flow Cyt, IHC-P	Human
242402	Anti-Cytokeratin 4 antibody [SP210] - BSA and Azide free	IHC-P	Human
242404	Anti-Slit2 antibody [SP122] - BSA and Azide free	IHC-P	Human,Mouse
242416	Anti-CTAG1B antibody [SP349] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
242417	Anti-EGFR (phospho Y1068) antibody [SP353] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
242418	Anti-ERK1 + ERK2 (phospho T202 + T204) antibody [SP327] - BSA and Azide free	Flow Cyt, IHC-P	Human
242419	Anti-eIF4G1 (phospho S1108) antibody [SP351] - BSA and Azide free	IHC-P	Human
242420	Anti-Lysozyme antibody [SP329] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
242422	Anti-Inhibin alpha antibody [SP323] - BSA and Azide free	IHC-P	Human
242423	Anti-Wilms Tumor Protein antibody [SP320] - BSA and Azide free	Flow Cyt, IHC-P	Human
242424	Anti-beta Catenin antibody [SP328] - BSA and Azide free	Flow Cyt, IHC-P	Human
242425	Anti-Wilms Tumor Protein antibody [SP321] - BSA and Azide free	Flow Cyt, IHC-P	Human
242426	Anti-Tcl1 antibody [SP248] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
242427	Anti-MUM1 antibody [SP114] - BSA and Azide free	Flow Cyt, IHC-P	Human
242428	Anti-TTF1 antibody [SP141] - BSA and Azide free	Flow Cyt, IHC-P	Human
242429	Anti-PAX8 antibody [SP348] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
242430	Anti-Lysozyme antibody [SP350] - BSA and Azide free	Flow Cyt, IHC-P	Human

243909	Anti-SOX2 antibody [SP76] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
243910	Anti-Carbonic Anhydrase 9/CA9 antibody [SP106] - BSA and Azide free	Flow Cyt, IHC-P	Human
243912	Anti-VEGF Receptor 2 antibody [SP123] - BSA and Azide free	IHC-P	Human
243914	Anti-PI 3 Kinase catalytic subunit alpha/PIK3CA antibody [SP139] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
243915	Anti-Glucose Transporter GLUT1 antibody [SP168] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
243920	Anti-Bcl6 antibody [SP155] - BSA and Azide free	Flow Cyt, IHC-P	Human
243921	Anti-MUC4 antibody [SP241] - BSA and Azide free	IHC-P	Human
243922	Anti-Cytokeratin 5 antibody [SP178] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
243923	Anti-Neutrophil Elastase antibody [SP203] - BSA and Azide free	Flow Cyt, IHC-P	Human
243924	Anti-LYPD3 antibody [SP204] - BSA and Azide free	IHC-P	Human
243925	Anti-CD16 antibody [SP175] - BSA and Azide free	Flow Cyt, IHC-P	Human
243926	Anti-CD13 antibody [SP182] - BSA and Azide free	Flow Cyt, IHC-P	Human
243929	Anti-Smad4 antibody [SP306] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human,Mouse
243930	Anti-Met (c-Met) antibody [SP44] - BSA and Azide free	Flow Cyt, IHC-P	Human
243931	Anti-Desmin antibody [SP138] - BSA and Azide free	ICC/IF, IHC-P, WB	Human,Mouse
243932	Anti-CD31 antibody [SP164] - BSA and Azide free	Flow Cyt, IHC-P	Human
243933	Anti-NRAS (mutated Q61 R) antibody [SP174] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
243934	Anti-CD10 antibody [SP179] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human,Mouse,Rat
243935	Anti-MAGEA1 antibody [SP188] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
243936	Anti-CD16 antibody [SP189] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
243937	Anti-PD1 antibody [SP269] - BSA and Azide free	Flow Cyt, IHC-P	Human
243938	Anti-SLC39A6/ZIP-6 antibody [SP300] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
243939	Anti-Integrin beta 3 antibody [SP324] - BSA and Azide free	IHC-P, WB	Human
243940	Anti-PTEN antibody [SP218] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
245729	Anti-Cyclin E1 antibody [SP146] - BSA and Azide free	IHC-P	Human
245730	Anti-TdT antibody [SP150] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
245731	Anti-CD3 antibody [SP162] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse
245732	Anti-CD23 antibody [SP163] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
245733	Anti-SPARC antibody [SP205] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
245734	Anti-Cyclin D3/CCND3 antibody [SP207] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
245735	Anti-Sialoadhesin/CD169 antibody [SP213] - BSA and Azide free	Flow Cyt, IHC-P	Human
245736	Anti-HLA-DP antibody [SP228] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
245737	Anti-Indoleamine 2, 3-dioxygenase antibody [SP260] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
245738	Anti-KMT6 / EZH2 antibody [SP129] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
245739	Anti-FOXP1 antibody [SP133] - BSA and Azide free	Flow Cyt, IHC-Fr, IHC-P, WB	Human
245740	Anti-p53 antibody [SP161] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
245741	Anti-CD13 antibody [SP187] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
245742	Anti-Sialoadhesin/CD169 antibody [SP216] - BSA and Azide free	Flow Cyt, IHC-P	Human
245743	Anti-VPAC1 antibody [SP234] - BSA and Azide free	Flow Cyt, IHC-P	Human
245744	Anti-Nucleophosmin antibody [SP236] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human,Mouse,Rat
245745	Anti-CD79b antibody [SP240] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
245746	Anti-CD103 antibody [SP301] - BSA and Azide free	Flow Cyt, IHC-P	Human
245747	Anti-PSMA antibody [SP29] - BSA and Azide free	IHC-P	Human
245748	Anti-Ku80 antibody [SP130] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
245749	Anti-Villin antibody [SP145] - BSA and Azide free	Flow Cyt, IHC-P, IP, WB	Human,Mouse
245750	Anti-Lambda Light chain antibody [SP147] - BSA and Azide free	IHC-P, WB	Human
245751	Anti-GIV (phospho Y1765) antibody [SP158] - BSA and Azide free	Flow Cyt, IHC-P	Human
245752	Anti-MAP4K4/NIK antibody [SP177] - BSA and Azide free	Flow Cyt, IHC-P	Human
245753	Anti-Calponin 1 antibody [SP180] - BSA and Azide free	IHC-P	Human
245754	Anti-Calponin 1 antibody [SP181] - BSA and Azide free	IHC-P, WB	Human
245755	Anti-MAGEA1 antibody [SP194] - BSA and Azide free	Flow Cyt, IHC-P, WB	Human
245756	Anti-CD134 / OX40L receptor antibody [SP195] - BSA and Azide free	IHC-P	Human
245757	Anti-Factor XIIIa antibody [SP196] - BSA and Azide free	IHC-P	Human
245758	Anti-C19orf2 antibody [SP215] - BSA and Azide free	Flow Cyt, IHC-P	Human,Rat
245759	Anti-Sall4 antibody [SP289] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
245760	Anti-SOX10 antibody [SP267] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P, WB	Human
246319	Anti-BFL-1/GRS antibody [SP435]	IHC-P	Human
246348	Anti-Calponin 1 antibody [SP193] - BSA and Azide free	IHC-P	Human
246349	Anti-LYPD3 antibody [SP208] - BSA and Azide free	IHC-P	Human
246350	Anti-PTEN antibody [SP227] - BSA and Azide free	IHC-P, WB	Human
246351	Anti-HLA-DPB1 antibody [SP229] - BSA and Azide free	Flow Cyt, ICC/IF, IHC-P	Human
246352	Anti-VPAC2 antibody [SP235] - BSA and Azide free	Flow Cyt, IHC-P	Human
246353	Anti-Fbxw7 antibody [SP237] - BSA and Azide free	Flow Cyt, IHC-P	Human
246709	Anti-CD42b antibody [SP219] - Low endotoxin, Azide free	Flow Cyt, IHC-Fr, IHC-P	Human,Mouse
246794	Anti-CD276 antibody [SP206] - Low endotoxin, Azide free	Flow Cyt, IHC-P, WB	Human
246798	Anti-CD16 antibody [SP189] - Low endotoxin, Azide free	Flow Cyt, IHC-P, WB	Human
246799	Anti-CD3 antibody [SP162] - Low endotoxin, Azide free	Flow Cyt, IHC-P, WB	Human
253247	Anti-BFL-1/GRS antibody [SP435] - BSA and Azide free	IHC-P	Human
15556	Anti-Phosphoserine/threonine/tyrosine antibody [SPM101]	ELISA, Flow Cyt, FuncS, ICC/IF, IHC - Wmt, IHC-Fr, IHC-P, IP, WB	Species independent
40084	Anti-Glucose Transporter GLUT1 antibody [SPM498]	Flow Cyt, ICC/IF, IHC-FoFr, IHC-Fr, IHC-P, WB	Mouse, Rat, Human
228461	Anti-BRAF (mutated V600 E) antibody [VE1]	IHC-P	Human (predicted: Mouse, Chicken)